

PC Login

Admin

password: 3000hanover

****Exporting chromatograms to CSV (for replotting on prism) ****

Follow this helpfull[link](#).

Or use the method that has been set up (2023-11-01): 3D...**CSV**.pmx

1. Open the data file you want to export
2. You can just opt to click the usual 3D UV....
3. In the **Home** tab, Under **Methods** → **Open Method** ** -> **3D...CSV.pmx** (**This methods saves the file in C:\CDSProjects\Installation\CSVexport\.... It will make a folder for each of your data file**) - **Make sure the processing method at the bottom left of the screen is set to this method. Click **Save. -> Reprocess All.**
4. Open the folder in C:\CDSProjects\Installation\CSVexport\ and transfer you file for plotting

From:

<https://doku.aephir.net/> - **SafaviLab Wiki**

Permanent link:

https://doku.aephir.net/doku.php?id=documents:agilent_hplc

Last update: **2023/11/01 09:25**

